

C8⁺: QUALITY PLUS FLEXIBILITY.

NEW
NOW ALSO AVAILABLE
WITH LOUNGE
UPHOLSTERY

SIRONA.COM

The Dental Company

sirona.

THE REAL FIT FOR YOU.

The Sirona C8⁺ dental treatment center offers unrivalled ergonomics, flexibility and esthetics, plus a wide range of options that enable you to customize it to your specific requirements. **Enjoy every day. With Sirona.**

ABSOLUTE RELIABILITY

Outstanding quality made in Germany,
30,000 installations worldwide

HIGH FLEXIBILITY

Three versions, numerous
individual configuration options

INTUITIVE OPERATION

Easily understandable
operating concept



ADAPTIVE TO CHANGING REQUIREMENTS.

With its three versions and the flexible configuration options, the C8+ cannot only be perfectly adapted to your individual treatment method but is also able to meet any changing requirements.



The whip arm version
The flexible spring-loaded joints of the whip arms increase the operating radius of the instrument hoses. The dentist element can thus be brought into an ergonomically perfect position and the molar areas are always easily accessible.



The OTP version
The variably positionable support arm can be adjusted in height independent of the chair position as it is fastened on the floor-mounted water unit. The dentist element provides lots of storage space.





The cart version
The cart is directly connected to the patient chair and is freely positionable. This version is just as easy to operate as the whip arm and OTP versions.



The turn option
Due to the symmetric design of the treatment center the C8+ can easily and quickly be converted from a right-handed system to a left-handed system.



The conversion option
And if you wish to sell your treatment center at a later date you can always very easily convert it from an OTP unit to a whip arm unit or vice versa.



MORE THAN ERGONOMICS.

Easier workflow, perfect results – with this goal, Sirona continuously develops its ergonomics program.

Sirona's ergonomics program is based on decades of experience in ergonomics. And as an innovator in linking digital technology with system solutions, we always think ahead. The ergonomics circle stands for Sirona's comprehensive perspective of the relationship of ergonomics and workflow in the daily practice routine. The round shape also indicates that Sirona's ergonomics program is being developed continuously.



TREATMENT WITH COMFORT.



The user comfort of the C8⁺ has its roots in the Sirona Ergonomics Program and can be further increased through options such as the four-way foot switch and the HUGO working stool.



Intuitive operation

Everything can be controlled easily and precisely via the operating panel. The memory function saves the last settings used.



Innovative instruments

Sirona instruments are widely known for their top performance in all categories and with regard to all treatment applications.



Handsfree operation

Either with the four-way foot switch or with the C⁺ foot control C8⁺ can be operated even more intuitively.



Optimum visibility

LEDview with three swivel axes for optimal positioning. Instruments also available with innovative LED light upon request.



INTUITIVE SITTING POSITION

More freedom of movement

Flat backrest for greater legroom

Optimized access

Perfect coordination and flexible positioning of dentist and assistant elements

Smart working stools

HUGO, CARL, PAUL



COMFORTABLE PATIENT POSITIONING

Perfect adjustment of the treatment center

Flexible height adjustment, ErgoMotion: synchronized movement of seat and backrest

Ideal support for the spine

No "stretching effect", Lounge upholstery



OPTIMUM VISIBILITY

Ideal illumination

LEDview, LED instrument light, SIROLUX E

Flexible headrests

Flat headrest, articulated headrest, MultiMotion headrest

Perfect intraoral images

SIROCAM 3 with 3 focus ranges



INTEGRATED WORKFLOW

Intuitive operation

Simply structured operating panel, C⁺ foot switch for handsfree operation of the chair programs

Modern equipment

Integrable camera, instruments with LED light

PERFECT RESULTS

Clinically perfect

Esthetically perfect

Perfect for the success of your practice



MultiMotion headrest

Can be inclined, tilted and rotated in virtually any chosen direction and can easily be operated with just one hand.



HUGO

The smart working stool which ensures you to adopt correct posture intuitively.



ErgoMotion

It is not necessary to readjust the headrest extension when reclining or inclining the backrest.

RELAXED PATIENTS.



C8+ not only offers a high level of comfort and convenience for you as a dentist, but also for your patients.

With its anatomically shaped upholstery, perfect support of the spine through ErgoMotion and the optional Lounge upholstery C8+ provides exceptional patient comfort. By means of the SIVISION 3 patient communication system you can rapidly gain even more trust from your patients. And the customizable choice of colors additionally creates an atmosphere of well-being.



NEW: Lounge upholstery
The soft surface increases comfort, shows style and allows for a relaxed treatment atmosphere merely through its appearance and texture. Available in the following colors: Carbon, Pacific and Rosewood.



SIVISION 3
Patient communication with intra-oral camera and 22" monitor as well as an optional direct connection to the practice network.

ErgoMotion
The combined motion of the backrest and the seat prevents stretching and compression of the spine.

BACKREST
The anatomic design of the backrest and headrest gives your patients a feeling of safety.

MODERN DESIGN

Sirona treatment centers have won numerous design awards. With its three color ambiances they blend in perfectly with every practice interior and create an atmosphere where patients feel comfortable.

THE ELEGANT WORLD

Classic, timeless combinations for an especially high-class look.



THE NATURAL WORLD

Soft, warm shades for a soothing treatment atmosphere.




THE VITAL WORLD

Vivid colors and bright tones for a modern, lively atmosphere.



PAINTWORK COLORS (CAN BE INDIVIDUALLY COMBINED)




With the color configurator at sirona.com you can see for yourself the effect of the colors you choose in your treatment center.

PERFECT HYGIENE.

The hygiene functions of the C8⁺ ease the burden on your team and make it easy to comply with all applicable infection control regulations.

The C8⁺ turns hygiene into an easy, safe and time-saving matter. This is not only ensured by the design of the treatment center with its smooth surfaces and removable parts but also by the sophisticated hygiene functions, such as, for example, the integrated anti-retraction valve. The integrated disinfection system with permanent disinfection and sanitization function provides for an excellent water quality.



Integrated disinfection system
For a consistently high water quality. Also helps to ensure DVGW approval and provides the back-flow preventer.



Removable parts
Hygiene-critical parts can be detached and thermo-disinfected or even sterilized.



Antistatic bellows
Dust-repellent and easy to clean, even with the chair being raised to its uppermost position.



Smooth surfaces
Especially easy to clean and thus extremely hygienic: the surfaces of the complete treatment center.



Touch pad cover
The silicone cover of the touch pad provides for an even higher hygiene level without impairing operation.

YOUR INDIVIDUAL OPTIONS.

| Options |
|--|
| Patient chair |
| Flat headrest |
| Articulating headrest |
| MultiMotion headrest |
| Standard upholstery |
| Lounge upholstery |
| Wide or narrow backrest on choice |
| Swiveling arm rest, left or right |
| 4-way foot switch for easy chair control |
| Valve set for external devices (incl. hydrocolloide support) |
| Booster cushion for small patients |
| Standard pneumatic foot control |
| C* foot control for dynamic micro motor control |
| Turn version with chair swivel for easy conversion for left/right handed treatment |
| Dentist element options (5 handpiece positions, removable instrument holder) |
| OTP support arm with hanging hoses (Traditional style) |
| Whip arm version (Continental style) 2 whip arm types on choice |
| Cart version |
| Standard 3-way syringe |
| 3-way syringe SPRAYVIT |
| Drive block for 2 drill drives |
| 3rd drill drive kit for additional turbine |
| Electric motor SL (halogen or LED light) |
| Electric motor SL ISO (halogen or LED light) |
| 6-hole turbine hose (halogen or LED light) for fiber optic air-driven instruments |
| Turbine hose Borden style |
| Piezo-ceramic scaler SIROSONIC L (halogen or LED light) |
| Intra-oral X-ray viewer |
| Satelec mini L.E.D. curing light |
| Intra-oral light source |
| Swiveling tray holder (for 1 or 2 standard trays on choice, only for whip arm version) |
| Silicone pad, provides non-slip work surface |
| Removable silicone sleeves for handle and touch pad |

| Options |
|--|
| Assistant element |
| Satelec mini L.E.D. curing light |
| Standard 3-way syringe |
| 3-way syringe SPRAYVIT |
| HVE |
| Saliva ejector |
| Touch pad for chair and cuspidor functions |
| Long support arm for two-handed treatment |
| Water unit options |
| Standard cuspidor bowl |
| Extended cuspidor for improved accessibility |
| Cover instead of cuspidor |
| Disinfection device for continuous waterline disinfection and sanitisation |
| Fresh water supply - bottle contents approx. 1.5 ltr – with option to switch over to public water supply |
| Water heater for instrument spray and cup fill water |
| Purge function – Automatic rinsing of instrument water lines |
| Preparation for wet suction with or without micro switches |
| Multi-station valve for wet suction systems |
| Dürr cuspidor valve for wet suction systems with centralized amalgam separator |
| Air-water separator for dry suction systems |
| Integrated amalgam separator according to ISO 11143 |
| SIVISION 3 |
| Preparation of dentist element support arm for camera |
| SIVISION 3 for camera in the dentist element and monitor on the light support column |
| Monitor 22" incl. speakers |
| SIROCAM 3 intra-oral camera |
| Cable set for PC connection |
| Stools |
| HUGO stool with optional foot release height adjustment |
| CARL stool with 360° release mechanism and backrest rotation |
| PAUL stool with 90° backrest rotation |
| Foot ring for PAUL stool |
| Working height for 1.50 – 1.80 m or 1.70 – 2.05 m of body size on choice |
| Operating light |
| SIROLUX E halogen light (22,000 Lux) |
| LEDview light with sensor switch |

YOUR EQUIPMENT OPTIONS.

In the following you find an overview of all available options regarding the equipment of your C8+ treatment center.



HUGO Freehand



HUGO Manual Plus



HUGO Manual



CARL Manual



PAUL stool



Operating panel for assistant



Swiveling armrest



Wide/narrow backrest



Booster cushion for small patients



Satelec mini L.E.D. curing light



Valve set for external devices



Fresh water bottle



SIVISION monitor 22"



SIROCAM 3



Intra-oral X-ray viewer



SIROLUX E

ALWAYS AT THE FOREFRONT OF INNOVATION!

As global innovation leader for dental equipment, we continuously invest in research and thus in the future of modern dentistry. By networking digital technologies with integrated solutions and optimizing the treatment workflow, we create improved treatment results, more comfort and safety for the patient as well as time and cost savings in everyday work. The combination of constant innovative power and globally growing sales and service structures makes Sirona the global market leader trusted by thousands of practices and labs around the world. **Enjoy every day. With Sirona.**



CAD/CAM systems

From pioneer to new standard. For 30 years we have been developing digital dentistry and creating new possibilities for the future practice and lab.



Imaging systems

Best image quality with the lowest dose. More than 100 years of developing x-rays for the dental practice make us the number 1 innovation partner.



Treatment centers

The business card of modern practices. We are striving to create the ideal ergonomic and innovative center. Individually tailored to the well-being and demands of the patient and dentist.



Instruments

Advantages that speak for themselves. We make sure that we provide the right balance of proven quality, individual ergonomics and innovative technology for user-friendly work.



Hygiene systems

Competence that gives you safety. When it comes to hygiene in the practice, we do not take any shortcuts.