

# CAN RELY ON.

All equipment options can be found at [SIRONA.COM/INTEGO](http://SIRONA.COM/INTEGO)



### INTEGO pro

- Integrated, comprehensive hygiene concept for optimal results.
- EasyTouch touchscreen and integrated foot control concept for easy, intuitive operation.
- Flexible assistant element suitable for two- or four-handed treatment.
- Numerous equipment options with a great value.
- Available as the TS model with hanging hoses or CS with whip arms.

# QUALITY YOU CAN RELY ON.

We don't just want to satisfy our customers, we want to impress them. This is why quality, performance, and service are key for us. Sirona treatment units combine high quality with real innovations that perceptibly enhance treatment comfort and workflow.



**Reliability and durability**  
Precise, elaborate quality tests and robust materials make INTEGO a reliable treatment center with an above-average service life.



**High quality**  
INTEGO immediately shows our attention to detail. The materials used – robust, shapeable plastics, solid steel and aluminum components on the inside, high-quality glass and upholstery – meet the high quality demands of a premium product.



**High standards**  
Sirona is a technology leader with more than 130 years of innovation history and sets high quality standards. Some 290 engineers and scientists in Research & Development not only pay close attention to quality; they also drive the innovation process.



**Future-proof design**  
At the Sirona innovation center, future trends are anticipated and products are prepared to be equipped for the technology demands of the future.



Subject to technical changes and errors in the text. Order No. A91100-M4-8812-01-7600, Printed in Germany, Dispo-No. 0 4601, DEW15 DD 0315.V1

# INTEGO. PERFORMANCE YOU CAN RELY ON.



SIRONA.COM

# GREAT VALUE YOU

Your work day is complex and places many different demands on the treatment center you choose. INTEGO gives you the perfect solution: Choose from two models that you can configure flexibly. Whether you select an INTEGO or INTEGO pro, a whip arm or a hanging hoses design – every option offers optimal performance at a fair price.



### INTEGO

- Treatment center with a compact water unit for optimal use of space.
- EasyPad user interface for easy, intuitive operation.
- Compact assistant element with good freedom of movement.
- Controlled water quality with a self-contained fresh water supply bottle.
- Available as the TS model with hanging hoses or CS with whip arms.

The Dental Company

sirona.

The Dental Company

sirona.

# INTEGO

INTEGO gives you the ideal workflow for your daily processes. is an extremely reliable partner for your practice.



**EasyPad**  
Clear and intuitive user interface.



**Fresh water bottle**  
For controlled water quality. Can also be filled with DENTOSEPT solution for permanent disinfection and manual sanitization.



**Automatic suction hose cleaning**  
Cleaning adaptor for automatic purging of saliva ejector and spray aspirator.



**Integrated torque control**  
The motor stops automatically when the set maximum torque is reached. This ensures greater safety and a perfect workflow for endodontics treatment.



**Assistant Element: Compact**  
The assistant element is very compact and ensures maximum freedom of movement. With the Compact model, you can work comfortably with an assistant.

# ERGONOMIC DESIGN YOU CAN RELY ON.




Every treatment is a complex procedure. Our ergonomic products provide optimal access to the patient and smooth treatment workflows with a high level of patient comfort. INTEGO combines every dimension of ergonomics to make treatment easier for you. Achieve perfect results even at the end of a long working day – with your expertise and the support of our treatment center. More about our ergonomic designed features at: [SIRONA.COM/INTEGO](http://SIRONA.COM/INTEGO)

## THE COLOR COMBINATIONS

### THE ELEGANT WORLD

Classic, timeless combinations for an especially high-class look.



### THE NATURAL WORLD

Soft, warm shades for a soothing treatment atmosphere.



### THE VITAL WORLD

Vivid colors and bright tones for a modern, lively atmosphere.



Scan the QR code and find out even more!  
[SIRONA.COM/INTEGO](http://SIRONA.COM/INTEGO)

# INTEGO pro

INTEGO pro gives you even more choice: All your requirements are covered thanks to the additional innovative operating and hygiene equipment options.



**EasyTouch**  
More operating comfort thanks to intuitive touchscreen technology and more operating and display options.



**ApexLocator**  
Can be integrated into the EasyTouch to allow permanent control of the working depth during endodontic treatment.



**Integrated cleaning adapter**  
For automatic purging of instrument hoses (AutoPurge) and sanitization of the water lines.



**Integrated disinfection system**  
Permanent disinfection of all water lines with DENTOSEPT. Automatic sanitizing function for regular intensive disinfection.



**Assistant Element: Comfort**  
The assistant element is very compact and ensure maximum freedom of movement. You can work comfortably even without an assistant.